

## Lab Notes

10/14 Lab Metabolic Tests (Ex 34-39)

### **Phenol Red (PR) Broth** -Durham Tubes **pg. 267**

- A simple broth that has a sugar in it (can vary)
  - PR **glucose, lactose & sucrose** (commonly used in lab)
- PR is a pH indicator
  - $\text{pH} < 6.8 \rightarrow$  yellow
  - $6.9 < \text{pH} < 7.3 \rightarrow$  red
  - $\text{pH} > 7.4 \rightarrow$  pink/magenta
- Carbohydrate fermentation results in lower pH (acidic end products) or acid end products + gas
- Gas will only be an end product if there is also acid

Results/Observations: turbidity  $\rightarrow$  growth, color  $\rightarrow$  pH, gas  $\rightarrow$  fermentation  
 $\rightarrow$  EXAM: understand what each result means & explain

### **Methyl Red (MR), Voges-Proskauer (VP)** **pg. 263**

- Methyl red test
  - 15 drops methyl red to MRVP broth tube
  - Mix
  - Instant color change
  - Methyl red is a pH indicator

Glucose fermentation  $\rightarrow$  stable acids

Results/Observations: Acid = red; (+) result

Basic = yellow; (-) result

no acid = no color

*S. epidermidis* (+)

- VP test (Do this test before MR Test!!!)
  - Transfer 1mL to another tube
  - Add 9 drops VPA (use blue nutcracker to break glass/release contents)
  - Mix (shake contents of tube)
  - Add 3 drops of VPB (use blue nutcracker to break glass/release contents)
  - Mix
  - Wait 15-20 min
  - Read every 10 min for 45 min
  - Do not shake tubes

Glucose fermentation  $\rightarrow$  acetoin + butanediol

After VPA and VPB

Results/Observations: red = (+) result; don't confuse with copper color

## **S. epidermidis (-)**

### **Nitrate test (NO<sub>3</sub>) -pg 271**

- Anaerobic respiration involving reduction of an inorganic molecule other than oxygen
  - N<sub>2</sub>O (nitrous oxide)
  - NO (nitric oxide) -gas, hormone, chemical signal: increases blood flow
  - NOS
- Nitrate reduction broth
  - Denitrification (converting nitrate to N<sub>2</sub>)
  - Nitrate → nitrite → nitric oxide → nitrous oxide → molecular nitrogen
    - This occurs when O<sub>2</sub> is not the final electron acceptor

#### **Nitrate test**

- Add 6 drops of nitrate reagent A and B
- Within 2 min, change in color occurs

Results/Observations: if red = (+) for nitrite (if nitrite is detected) NO<sub>3</sub> was reduced to NO<sub>2</sub> (or nitrate was reduced to nitrite meaning denitrification took place) TEST done if red!

If no red color, add a pinch of zinc powder

- Within 2 min, color change occurs
- Zinc + NO<sub>3</sub> = Red (nitrogen reduction did not take place) = (-) test result for NO<sub>3</sub> reduction
- If no red after zinc, (+) test because nitrate went through nitrogen reduction to a product further than nitrite

### **Citrate Test pg. 269 -slant (stab and smear)**

- Used to differentiate organisms based on ability to grow when essential nutrient (like carbon) is limited
- Sodium citrate = sole C source
- Ammonium phosphate = sole N source (nitrogen used for proteins/amino acids, nucleic acid -DNA; nucleotides; nitrogenous base)
- Citrate → alkaline (NH<sub>3</sub>) + ammonia products + Bromthymol blue (blue when pH inc/green when pH is neutral)
  - Ammonium is basic (causes pH to go up)

Results/Observations: blue when pH increases (BASIC)

Green when pH is neutral

S. epidermidis (-) -appeared green with little to no growth on slant

### **Tryptic Soy Agar (TSA) -slant**

#### **Toxic buildup**

- The ETC can lead to a buildup of Hydrogen Peroxide (H<sub>2</sub>O<sub>2</sub>)
  - 2H<sub>2</sub>O<sub>2</sub> ---- (via catalase) ----> 2H<sub>2</sub>O + O<sub>2</sub>

- In order to engage in cellular respiration, bacteria need to make/have catalase
- Catalase test
  - Use glass slide
  - Use applicator stick to smear your cells
  - Add 1-2 drops of H<sub>2</sub>O<sub>2</sub>
  - Look for bubbling (determines if catalase is present or not)
  - *S. epidermidis* bubbled immediately (+) test
- Oxidase test
  - Tests for cytochrome oxidase -enzyme in ETC
    - Take sterile cotton swab to pick up cells
    - Apply 2 drops of oxidase reagent
    - Within 2 min, color change occurs

Purple color = (+) result (cytochrome oxidase is present)

No color change (-) result (cytochrome oxidase is not present)

*S. epidermidis* = (-) test no color change occurred; cytochrome oxidase is not present

- **Spirit Blue Test(Fat Hydrolysis) pg. 276**
  - Lipase (exoenzyme)
    - Break down lipids
    - Lipids (too big to enter cell) —(lipase) —> C,H (smaller molecules that can enter cell)
    - Look for halos around colonies (transparent blue)
    - Use stock (slant) to inoculate spirit blue plate with just a straight line
    - *S. epidermidis* shows a halo around the colony; conclusion *S. Epidermidis* appeared to be capable of making the enzyme lipase; typically it shouldn't but somehow it does show a (+) test
    - No you got a negative test which was supposed to happen
- **Starch Hydrolysis test pg. 275**
  - Starch is a polysaccharide made of glucose monomers
  - Some cells can breakdown starch to glucose (requires an enzyme)
    - Amylase (some bacterial cells make it and some don't)
  - Starch is too big to enter cell, needs to be broken down via amylase into glucose
  - Starch —(amylase)—> glucose
  - Inoculate starch hydrolysis plate with a straight line
  - Add iodine -dark blue/black in presence of starch.
  - After addition of iodine; No zone of hydrolysis around streak; (-) test; does not create amylase

- **Casein Hydrolysis Test p. 275**

- Casein -protein found in milk
- Can be broken down by casease
- Casein — (casease) —> peptides/nitrogen behaving compounds (small enough to get into cell)
- Bacteria get smaller peptides; source of amino acids
- Inoculate starch hydrolysis plate with a straight line
- No white appearance of white halo around bacteria; *S. epidermidis* (-) test result

- **Mannitol Salt Agar pg. 387**

- Ingredients include: beef extract, peptone, NaCL (75g in 1L H<sub>2</sub>O = 7.5%)
  - Typical [salt] in bacterial cell = 0.85%
- Mannitol (fermentable; get acids)
- Phenol red (pH indicator)
  - Acidic = yellow
  - Neutral = red
  - Basic = hot pink
- Which ingredient is selective? -salt (not a lot of bacterial cells will be able to tolerate the salt concentration)
- Which ingredient is differential? -mannitol (not an essential ingredient, but allows us to differentiate the cells via color; doesn't restrict in growth or anything)
- A differential and selective medium
- Favors growth of staphylococcus
  - Differentiates: pathogenic staphylococcus (ferments mannitol) from nonpathogenic staphylococcus (does not ferment mannitol)
- Streak the slant
- Appears bright pink and orange reddish towards the bottom; small amount of growth = able to tolerate salt concentration → basic; more alkaline; tolerant of salt = (-) MSA test

- **EMB Agar (pg. 414)**

- Used to detect coliforms (fecal contamination)
  - gram negative, bacillus (rod shaped), lactose fermentation
- Selective medium -encourages growth of gram (-) organisms
  - Inhibits growth of gram (+) organisms
- Ingredients include: lactose, color indicators: sodium sulfate and basic fuchsin
  - Detects changes in pH

- Also double as gram (+) inhibitors
  - Inoculate plate with wavy streak
  - Little to no growth on plate; colorless/color of medium; lactose nonfermenters
  - Red with gold metallic sheen -indicates vigorous lactose fermentation; indicative of coliform bacteria (fecal contamination)
  - *S. epidermidis* is Gram (+) → (-) test
- **Urease test pg. 276**
  - Urea is a nitrogen containing compound; broken down by urease
    - Products: carbon dioxide and NH<sub>3</sub> ammonia (contributes to an alkaline pH)
  - Urease is a hydrolytic exoenzyme
  - If rxn takes place, pH should increase = basic)
    - Basic = bright pink
    - Acidic = yellow
  - Inoculate culture broth tube; just drop cap on tube; do not push down
  - Organism appeared bright pink = basic = (+) test
- **Indole Test (Tryptophan Degradation) pg. 276**
  - Inoculate culture broth tube w/loop full of cells
  - Tryptophanase is enzyme needed to
  - Tryptophan (amino acid)—(tryptophanase) → indole + pyruvic acid
  - Add 15 drops of Kovac's Reagent
    - Yellow/no color change = (-) test result
    - Pink/red = indole (+) test result; makes enzyme tryptophanase
    - Indole is byproduct therefore if rxn occurred, then indole (pink/red color) should appear
    - *S. epidermidis* = (-) result; no change in color occurred
- **Phenylalanine Deaminase pg. 277**
  - Phenylalanine(amino acid) — (phenylalanine deaminase) → phenylpyruvic acid + ammonia (bi-product)
  - Streak the slant w/loop
  - Add 6-8 drops of Ferric chloride added to slant
    - Detects phenylpyruvic acid → Green color appears in 5-10 min
    - Yellow (remained the same color) = negative test = no phenylpyruvic acid is detected
    - *S. epidermidis* remained same color after addition of acid (-) test
- **Motility pg. 132 ex. 17**
  - Semisolid motility agar -allows cells to move if capable
  - Use needle to stab center of the media (do not go all the way down) using your organism
  - Motility test has TTC -Turns red when reduced

- Organisms shows some signs of movement; diffuses away from line (+) motility
- **Pg. 137 (use organisms provided to view under microscope)**
  - Loop full of cells on slide; use toothpick and vaseline; place vaseline on corners of cover slip; place cover slip using vaseline to stick
- **KIA (Kligler Iron Agar) pg. 279**
  - Glucose -1g in 1L
  - Lactose & Sucrose -10g each in 1L
  - Phenol red
  - Sodium thiosulfate
  - Ferrous sulfate
  - Inoculate slant by stabbing with needle & streaking slant
  - Glucose is preferred sugar; monosaccharide; easy to breakdown; but its only 0.1% -so it runs out within 12 hours
  - 0.1% glucose → acid (small amount); then goes to alternate source
  - Beef extract/peptone → NH<sub>3</sub> -ammonia (raises pH)
    - Reversion -slant reverts, but not the butt
    - Basic slant, acidic butt (red slant yellow butt)
    - Glucose only fermenter results: basic slant/acidic butt, H<sub>2</sub>S (+/-) -may or may not have, gas (+/-) -may or may not have
  - Glucose and lactose and/or sucrose fermenter
  - 1% lactose, 1% sucrose (10x) more than glucose
  - lactose/sucrose → pH drops -acid (large amounts)
    - No reversion -Tube is all yellow
    - acidic/acidic, H<sub>2</sub>S (+/-), gas (+/-)
  - Cysteine (amino acid) + H<sub>2</sub>O + H<sup>+</sup> — (cysteine desulfurase) —> pyruvate
    - Byproducts include H<sub>2</sub>S and NH<sub>2</sub>
    - Sugar fermentation produces acids (H<sup>+</sup>)
    - Cysteine desulfurase works only in an acidic environment
      - Breakdown of Cys to pyruvate
      - Produces hydrogen sulfide (H<sub>2</sub>S) gas -stinks
    - H<sub>2</sub>S + Fe → FeS (ferric sulfide), a black precipitate
    - Results: any black is (+) for H<sub>2</sub>S

Red = basic;

yellow; acidic;

presence of gas = bubbles; slant moved to the top;

presence of H<sub>2</sub>S = black precipitate; stinky → ONLY works in acidic environment

## Plaque Assays Ex. 21

- Bacteriophages -virus that specifically targets bacterial cells
- The Lytic Cycle
  - Virus attaches to cell wall
  - Injects viral genetic info into cell
  - Directs cell to make more virus
  - Cell lyses
  - Virus released
- Plaques -clearings in lawn of bacteria
  - Where cells have been infected and died
  - One plaque was the result of one virus infecting one cell
  - You can count the number of plaques to determine the original phage concentration
- Conc D.F. = DF(dilution factor) x Volume Plated
  - DF(dilution factor) = depends on plate; A =  $10^{-1} \times 10^{-1} = 10^{-2}$ 
    - A:  $10^{-2}$
    - B:  $10^{-3}$
    - C:  $10^{-4}$
  - Volume plated =  $10^{-1}$ 
    - Countable plaques range from 25-250
- PFU(plaque forming unit)/mL = (# plaques/FDF)

## Soil Antibiotics

Obtain dirt in bag

Tare paper; weigh out 1g of dirt

Add 1g dirt into 99ml bottle of water = 100 fold dilution of dirt

Mix & let settle

Use disposable pipette to draw up 1mL of liquid

Transfer to 1st plate & spread

Transfer 0.1mL to plate & spread (about 2 and a half drops)

Let plates sit right side up

Leave at room temperature

Results: should expect bacterial growth; more so on 1.0ml plate

## Gram Stain Unknown -Gram Negative Bacilli

Metabolic tests -flow chart

1. Citrate Test

11/6/25

Ex. 31: Soil Antibiotics (cont'd)

- Observations:
  - 1ml: lots of colonies; many different shapes/sizes; crescent shapes colonies
  - 0.1ml: significantly less than 1mL, as expected
- Next Step: Inoculate loop with a single colony and add one streak onto a new plate (1 per pair)

Ex: 56 Epidemic

- Candy #23 → inoculate bacteria and shake hands w 4 people; inoculate plate after each person
- 

11/18 LAB

Ex 131: soil Abx

Ex 43: Bacterial Food Counts

Ex 60: ELISA

Enzyme-Linked Immunosorbent Assay (ELISA)

- Direct (Sandwich) ELISA
  - Detects Ag in patient sample
  - Start with bound Ab
  - Add sample
  - Add enzyme linked Ab
  - Both Abs bind different epitopes
  - Add substrate
- Indirect ELISA
  - Detects Abs in patient sample
  - Start with bound Ag
  - Add sample
  - Add enzyme linked Ab
  - Add substrate

→ sample c & e appeared clear → negative result

Ex 61: Blood Typing

Blood Types

- Detects presence of A and/or B antigens
- Surface of RBC is covered with these proteins
- Helps with identification and signaling
- You develop antibodies against the antigens that your RBC lacks



#### Agglutination Test (test subject: 20)

- Blood typing
- Add anti-A, anti-b, anti Rh antiserum (separately) to blood
- Look for agglutination
- Pg 452

→ AB negative blood type

#### Ex 54: Blood Agar

##### Blood Agar Plates (BAP)

- Tryptic Soy Agar (TSA)
- 5% sheep's blood
- Swab your throat

#### Streptococci

##### Lancefield Classification system for streptococci

- Group A strep
  - Streptococcus pyogenes
  - Cause pharyngitis, scarlet fever
- Group B Strep
  - Streptococcus agalactiae
  - Part of normal flora of gut and female urogenital tract
- During delivery baby can acquire and develop neonatal sepsis, meningitis

##### NonLancefield Groups

- Viridans: Streptococcus mutans
  - Found in human oral cavity
  - Contributes to tooth decay
  - Adheres to tooth surface
  - Metabolizes sucrose to glucan dextran (sticky)
  - Biofilm forms (other cells collect)
  - Sugars fermented
  - Produce acid
  - Tooth decay

##### Differentiation of streptococci

- Group A and Group B
  - B-hemolytic
  - Complete destruction of RBC
  - Produces a clear zone on Blood Agar Plate
- Viridans and Strep. Pneumoniae
  - A-hemolytic
  - Incomplete destruction of RBC
  - Greenish discoloration

- Y-hemolytic
  - No apparent RBC damage
  - Enterococcus faecalis

#### Ex 15: Endospore Staining

- Endospore -cell formed within cells of certain bacteria
- Extremely resistant to adverse conditions
  - Temperature
  - Desiccation
  - UV light
  - Lack of nutrients
  - pH
- Can be formant for a long time
- Bacterial growth curve
  - Towards end of stationary phase, sporulation occurs

#### Endospore stain → differential stain

- Use Bacillus subtilis
- Make smear, air dry, and heat fix slide
- Flood with malachite Green (basic)
  - Use paper towel
  - Keep moist
- Steam **10 minutes**
  - Allows dye to penetrate spores
- Rinse thoroughly with water
- Counterstain: Safranin (basic) for **5 min**
- Water rinse
- Oil immersion

#### Acid Fast Stain Ex 16

- Tests for mycolic acids in the cell wall
- Mycolic acid -the waxy mycolic acid repels typical aqueous stains
- Carbofuchsin is a lipid soluble stain, and thus penetrates the waxy cell wall
- Staining enhanced by steaming

#### Procedure

- Prep bacterial heat fixed smear with Mycobacterium smegmatis AND E. coli
- Add carbofuchsin **12 min**
- Rinse
- Add acid-alcohol **15 seconds**
- Rinse
- Counterstain: Methylene blue **60 seconds**
- Rinse
- Oil Immersion

Cells that have mycolic acid will be magenta/reddish (acid-fast)

Cells with no mycolic acid will be blue (nonacid-fast)

E. Coli (nonacid-fast)

M. smegmatis (acid-fast)

### Capsules

- Protective polysaccharide layer
- Makes cell antiphagocytic
- Prevents desiccation

### Capsule Stain Ex 13

- Bacillus megaterium
- Prepare smear
- Air dry (NO HEAT FIX -you'll shrivel the capsule)
- 1% crystal Violet for 2 min
- Rinse with Copper Sulfate
- Oil Immersion

### Capsule Stain Ex 13 (mod)

- Bacillus megaterium
- Make negative stain using Congo Red (add congo red to edge, add cells, push/pull method)
- Air dry (NO HEAT FIX -you'll shrivel the capsule!)
- Counterstain: Maneval's Stain for 2 min
  - Background becomes dark blue and cells remain pink
  - Maneval's Stain is acidic
  - Turns Congo Red blue
  - Stains cells pink
- Rinse gently to rid excess Maneval's stain
- Oil immersion